---

## **ENGLISH**

The instructions for use of orthodontic appliance (EXALIGNER®) – custom-made device – produced by EXAKTUS are described in this document. The complete treatment is approved by a dentist or orthodontist responsible for the prescription (via a digital platform) and controlled by national laws and regulations.

### **1. INTENDED USE**

The EXALIGNER® is a medical orthodontic aligner made clear-transparent and abrasion-resistant copolyester material as well as thermoplastic elastomer. It is prescribed for malocclusion treatments and controlled by dentists and orthodontists. The aligners along time and with a correct and continuous use will apply some pressure in the teeth to adjust the position.

The aligners must be placed on the teeth and worn for at least 22 hours every day. In the beginning, it is normal to feel pressure applied to the teeth.

The patient should always follow the recommendations prescribed by the doctors and use the aligners in a responsible way.

### **2. MATERIALS INCLUDED**

Each box contains 1 full treatment per prescription (for 1 patient) according to the information approved by the professional team processing the initial information submitted via the digital platform.

### **3. INSTRUCTIONS**

- Place the aligners on top of your teeth and snap them on by applying balanced pressure on both sides of aligner.
- Aligners fit tightly and you may feel some pressure on your teeth, which is normal and is required to move teeth.
- After insertion of aligner, you may notice a slight “lisp” which will go away in few days, as your tongue adjusts to having aligners in your mouth.
- Aligners should be worn day and night for at least 22 hours a day, except while eating, brushing, and flossing.
- Avoid chewing gums and hot liquids while wearing aligners.
- You may feel “Dry mouth” at the start of the treatment. Drink plenty of water.
- Wear each aligner for two weeks or as instructed by your doctor.
- Most of the tightness will go away at the end of 1 week but if the aligners are still too tight, then do inform your doctor about it.
- Clean aligners everyday using toothbrush and normal tap water. Cleaning tablets can also be used to clean aligners. Do not use hot water.
- Avoid using toothpaste and mouthwash to clean aligners.
- When not in use always store aligners in the aligner box provided.
- In case of loss, damage, or breakage of your aligners, please report to your treating doctor immediately.

### **4. WARNINGS AND PRECAUTIONS**

- The aligners should be used in accordance with the treatment.
- In case of discomfort, fitting issues or other problems during use please contact your dentist or orthodontist.
- Remove always for meals and clean with a toothbrush and water.

### **5. POTENTIAL CONTRAINDICATIONS**

- Any type, method, or movement indicated as difficult in the case complexity guidelines.
- Centric-relation and centric-occlusion discrepancies.
- Teeth with short clinical crowns.
- Arches with multiple missing teeth.
- Patients whose second molars have not yet erupted.

- Patients with poor oral hygiene.
- Patients with active periodontal disease.
- Patients with dental prosthetics or implants.
- Patients with severe bruxism.

## 6. STORAGE AND STABILITY

- Store at room temperature, away from direct sunlight and moisture sources.
- Do not stack heavy boxes or objects directly over the box or over the aligners itself.
- Handle with care.
- Disposal must be made according to official regulations.

## 7. SYMBOLS

Symbols	Description	Symbols	Description
	Manufacturer		Batch code
	Consult instructions for use		Reference
	Medical device		Non-Sterile
	Caution		Date of manufacture
	Patient Nr	RX Only	Just with doctor's prescription

## 8. OTHER INFORMATION

The EXALIGNER® is produced in full compliance with the information submitted via the digital platform by the dentist and orthodontist responsible for the respective patient.

## 9. MANUFACTURER'S CONTACT INFORMATION

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In case of complaints or adverse incidents, you should contact the competent authority of your country and EXAKTUS.

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